




10 PL



Fractional IPL

—

- Blemish · Freckle · Age spot · Pigmentation
- Cafe au Lait · Acne · Partial lentiginosis · Flushing

FDA

CE

MFDS
KOREA

10PL is an IPL device that significantly lowers the risk of burns by applying the fractional treatment concept to the existing IPL. Also, the tip size is 45mm X 18mm, which allows of quick treatment of large areas.

10PL Application

*To be released





- 415nm** Melasma Epi, Acne
- 570nm** Melasma Dermi
- 590nm** Facial Veins, Rosacea




- 530nm** Vascularity, Pigmentation
- 700nm** Hair Removal
- 830nm** Skin Tightening

10PL Advantages

Conventional IPL



10PL




168 celles per irradiation area
(1.4mm x 1.59mm)

10PL Advantages


1) Treat Blemish using 35J High Power Fractional Beam

2) Select Wavelength using Turn Dial

Treat various lesions with three wavelengths without replacing filters!




3) 45mm x 18mm tip size reduces treatment time



4) Fast Procedure with Beam Irradiation Speed of 1Hz or less

Before & After



10 PL

Fractional IPL

This product is an aesthetic device of medical grade,
please refer to the precautions and treatment instructions before using it.



Tentech Inc. 4F, 5F Namjeon Building, 326, Bongeunsa-ro, Gangnam-gu, Seoul, Korea
TEL. 02-6954-2203 | FAX. 02-3444-4999 | www.tenlaser.com | For Export Only