
# Applicable to

- · Looking for a natural lifting
- Those who want skin care without interrupting daily life
- · Concerned about decreased skin elasticity
- · Wanting more meticulous skin care

# **Cartridge Types**

	Cartridge	Area
FACE	1.5mm	Eye rim, forehead
	2.0mm	Cheekbone, cheek
	3.0mm	Double chin, neck, cheek
	4.5mm	Cheek
BODY	7.0mm	Arm
	10.0mm	Thigh
	13.0mm	Abdominal





This product is not a medical device. It is a reference material to help understand the product.



**ONE THERA** 

DUAL MODE ULTRASOUND SYSTEM


### What is One Thera?

One Thera is a device that delivers concentrated thermal energy to the dermis and fat layer of the skin through various types of cartridges. It can be used for both FACE & BODY, and it is recommended to select a cartridge and mode suitable for the area to be used and operate in an optimal way.



## Long Pulse Mode?

Technology that divides one pulse into multiple pulses to deliver deep heat and further reduce pain.



### OTHERS

## Normal mode & Long pulse mode

Normal: Delivers thermal energy to the desired skin layer Long pulse mode: Transmits thermal energy deeper

- Approximately 4 minutes for 3,000 shots up to 12Hz speed
- Various applications possible by using two modes



#### 01 When can I see the results?

Normal: Delivers thermal energy to the desired skin layer Long pulse mode: Transmits thermal energy deeper

### 02 How often can I receive it?

Depending on the individual, you can expect immediate results immediately after use, and you will see better results after 2 to 4 weeks.

### 03 How is the pain level?

Depends on the individual, but generally it is a tingling degree of pain. If you are sensitive to pain, lower the energy level.

### 04 What are the precautions?

If the skin is red, it usually subsides on the day of treatment. It is recommended to refrain from sauna. drinking, and smoking for 3 to 5 days after the procedure.