

Regevan[®] fill
Fine



















Rhino[®] fill
Premium



TEN[®]TECH


Tentech Item

Regevan® fill Fine Rhino® fill Premium

									
<p>10 THERMA</p> <p>Monopolar RF</p>	<p>10 THERA</p> <p>2 Line HIFU</p>	<p>10 PL</p> <p>168 Cell Fractional IPL</p>	<p>10 CELL</p> <p>3cm x 3cm CO2 Fractional Laser</p>	<p>ONE THERA</p> <p>Dual Mode Ultrasound System</p>	<p>TEN HI</p> <p>Fractional & Thermal RF</p>	<p>10 SONO</p> <p>DSS & NBS Dual Mode System</p>	<p>LIPO THERA</p> <p>Dual Mode Ultrasound System</p>	<p>10 LUX</p> <p>High Intensity & Multiple Light Therapy Device</p>	<p>10 LUX HAIR</p> <p>PhotoBioModulation Device</p>
			 <p>Rhino fill Premium</p> <p>Regevan fill volume MAX</p> <p>Regevan fill volume PRO</p> <p>Regevan fill volume UP</p> <p>Regevan fill Fine</p> <p>Monophasic Filler [for Export]</p>			<p>Regevan[®] Pn</p> <p>Polynucleotide Sodium</p>			
<p>Regevan[®] LED Mask</p> <p>Dr. Oracle LED Mask 4,474 LED Lights</p>	<p>10 mono Thread</p> <p>Mono Thread [Registered for NMPA]</p>	<p>10 cog Thread</p> <p>Asymmetrical Cog Thread</p>					<p>Regevan[®] C</p> <p>Growth Factor Ointment</p>	<p>Regevan[®] S</p> <p>Growth Factor Spray</p>	<p>Regevan[®] ALO CARE</p> <p>Growth Factor Ampoule</p>

Regevan[®] fill Fine Rhino[®] fill Premium

For temporary improvement of facial wrinkles, it has three products:
Regevan Fill, and Rhino Fill
which were developed using world-class MDM technology,
and contains lidocaine to reduce pain during treatment.







Key technologies for manufacturing hyaluronic acid filler

Regevan® fill Fine Rhino® fill Premium

Class	Tentech Inc.	G*	A*
Country	Korea	Switzerland	Ireland
Brand	Regevan® fill <small>Fine</small> Rhino® fill <small>Premium</small>	Re*	Ju*
Core Tech.	MDM	NASHA	Hyla-cross
Phase	Mono Phasic	Bi Phasic	Mono Phasic
Remark	Recognized worldwide for its unique core technology	[How most filler companies are applied] · Cross-link : Short cross-link using BDDE	

Cross - link Comparison

Regevan® fill Fine Rhino® fill Premium

	Regevan® fill <small>Fine</small> Rhino® fill <small>Premium</small>	Normal Filler
Cross Linking Agent	DVS (Divinyl Sulfone)	BDDE (1,4-Butanediol diglycidyl ether)
Structure		
Mass (Dalton)	118.15	202.25
Molecule Length (Ratio)	1	3
Number of atoms forming a bridge between HA molecules	5	14 (Possibility of surplus atoms ↑)
Characteristics	Difficult process / Densely cross-linked	Simple process / Loosely cross-linked
Cross-link Completion	Completely cross-linked	'Pendant' phenomenon caused by incomplete crosslinking
Cross-link Illustration		

PENDANT Phenomenon in Filler Products Using BDDE as a Crosslinking agent

Determination of modification degree in BDDE-modified hyaluronic acid hydrogel by SEC/MS

Biao Yang^{*a}, Weuping Guo^{*a,b}, Hengchang Zang^{*a}, Jianjian Liu^{*b}
a School of Pharmaceutical Science, Shandong University, Jinan, China b Bloomage Freda Biopharm Co, Ltd, Jinan, China




Figure 2. When dissolved in water, hyaluronic acid (HA) behaves as a fluid, with excellent biocompatibility but poor mechanical properties (A). Modification of HA molecules by cross-linking improves mechanical properties by creating gels that have a firmer structure and are able to resist degradation (B). Modification does not necessarily cross-link HA to other HA molecules, resulting in a pendant cross-linker (C). Such structures often result in softer gels.

ORIGINAL ARTICLES

Comparative Physical Properties of Hyaluronic Acid Dermal Fillers



EFFERY KABLIK, *GARY D. MONTHIEIT, MD, LIPING YU, PHD, GRACE CHANG, AND JULIA GERSHKOVICH

Many BDDE technology imitations face Pendant challenges

Physical Property Composition Technology

Regevan® fill Fine Rhino® fill Premium

HA Filler Physical Properties Comparison (Monophasic vs Biphasic)

Type	MONO PHASIC	BI PHASIC
Illustration		
Physical Property	Type 1 HA (GEL)	Type 2 HA (GEL + LIQUID)
Product Characteristics	<ul style="list-style-type: none"> · High viscosity · Low elasticity 	<ul style="list-style-type: none"> · Low viscosity · High elasticity
Example	Regevan Fill, Rhino Fill, Fine Fill, Juvederm, Neuramis, Chaeum	Restylane, Yvoire, Cutegel

Comparison of Major Fillers


Regevan® fill Fine Rhino® fill Premium

Sortation	MDM HA Filler	General HA Filler	Chemical Filler
Persistence	0	X	0
MOLDING ability	0	X	0
Viscoelasticity	0	X	0
Safety	0	0	X
Biocompatibility	0	0	X
Biodegradability	0	0	X

MDM Filler™ is

Globally recognized for its safety and efficacy,
and it lasts longer with its excellent molding power.





MBMT™

Micro Bead Monophasic Technology
Integrated Microbead to manufacture Mono-phasic fillers

► Adjustment of viscoelasticity according to the purpose of use



DMCT™

Divinyl Sulfone Multi-Crosslink Technology
Fully cross-linked with Multi-cross-link using DVS

► Excellent viscoelasticity and molding power



MDAP™


Multiple Degree Amphiphilic Purification
Application of Multi-Step Refining Technology

► Manufacturing pure and secure products





MDM Filler™ Manufacturing Method

Regevan® fill Fine Rhino® fill Premium



Specification

Regevan® fill Fine Rhino® fill Premium

Product	Regevan® fill Fine	Rhino® fill Premium
Image		
Category	Medical Device	
HA Concentration	20mg/mL	
Viscosity	★★★★☆	★★★★★
Volume	1mL X 1 Syringe	
Lidocaine	0.3%	
Cross-link Completion	DVS	
Areas	All other areas	Nose

Regevan[®] fill Rhino[®] fill

Fine Premium

All information in this document is subject to copyrights and other intellectual property rights of Tentech Inc. and other licensees. No modification, reproduction or copy of any or any part of this material may be made without prior written permission from Tentech Inc. and licensees. Additional information is available upon request. The following are trademarks or registered trademarks of Tentech, Inc. in Korea and other countries.

